

PCR targeting system in *Streptomyces coelicolor* A3(2)

Bertolt Gust, Tobias Kieser and Keith Chater, John Innes Centre, Norwich Research Park, Colney, Norwich NR47UH, UK, Tel: +44 (0)1603 452751 Fax: +44 (0)1603 456844

Introduction

Many bacteria are not readily transformable with linear DNA because of the presence of the intracellular *recBCD* exonuclease that degrades linear DNA. However, the λ RED (*gam*, *bet*, *exo*) functions promote a greatly enhanced rate of recombination when using linear DNA. By exploiting this, Datsenko and Wanner (2000) made 40 different disruptions on the *E. coli* chromosome by replacing the wild-type sequences with a selectable marker generated by PCR using primers with 36 nt homology extensions.

The strategy for PCR-targeting for mutagenesis of *Streptomyces coelicolor* is to replace a chromosomal sequence within a *S. coelicolor* cosmid (Redenbach *et al.*, 1996) by a selectable marker that has been generated by PCR using primers with 39 nt homology extensions. The inclusion of *oriT* (RK2) in the disruption cassette allows conjugation to be used to introduce the PCR targeted cosmid DNA into *S. coelicolor*. Conjugation is much more efficient than transformation of protoplasts and it is readily applicable to many actinomycetes (Matsushima *et al.*, 1994). The potent methyl-specific restriction system of *S. coelicolor* is circumvented by passaging DNA through a methylation-deficient *E. coli* host such as ET12567 (MacNeil *et al.*, 1992). Vectors containing *oriT* (RK2; Pansegrau *et al.*, 1994) are mobilisable *in trans* in *E. coli* by the self-transmissible pUB307 (Bennett *et al.*, 1977, Flett *et al.*, 1997) or the non-transmissible pUZ8002, which lacks a *cis*-acting function for its own transfer (Kieser *et al.*, 2000).


To adapt the procedure of λ RED mediated recombination for *Streptomyces*, cassettes for gene disruptions were constructed that can be selected both in *E. coli* and in *Streptomyces* (Table 1). After a single disruption with an *oriT*-containing cassette, further disruptions can be performed on the same cosmid using *oriT*-free cassettes containing alternative selective markers. The λ RED recombination plasmid pKD20

(*E. coli* Genetic Stock Center CGSC Strain # 7637) was modified by replacing the ampicillin resistance gene *bla* with the chloramphenicol resistance gene *cat*, generating pIJ790, to permit selection in the presence of Supercos1-derived cosmids (ampicillin and kanamycin resistance).

Name of plasmid	Resistance-marker	Resistance	Concentration for <i>E. coli</i>	<i>oriT</i>	Size of template
pIJ773 Fig. 5	<i>aac(3)IV</i>	apramycin	50 µg/ml LB	+	1382 bp
pIJ778 Fig. 6	<i>aadA</i>	spectinomycin streptomycin	50 µg/ml LB 50 µg/m LB	+	1425 bp
pIJ779.	<i>aadA</i>	spectinomycin- streptomycin	50 µg/ml LB 50 µg/ml LB	-	1057 bp
pIJ780 Fig.7	<i>vph</i>	viomycin	30 µg/ml DNA	+	1497 bp
pIJ781	<i>vph</i>	viomycin	30 µg/ml DNA	-	1622 bp

Table 1: Disruption cassettes containing different resistance markers with and without *oriT*: All disruption cassettes were cloned into the *EcoRV* site of pBluescript SK II (+) allowing the isolation of a *EcoRI/HindIII* fragment for use as template for the PCR reaction. The size of the cassettes includes the 19 bp and 20 bp primer site (see section 2: “primer design”) which are identical in all disruption cassettes. The resistance genes with or without *oriT* are flanked by FRT sites (FLP recognition targets) which allows FLP-mediated excision of the cassette (see section 7: “FLP-mediated excision of the disruption cassette”).

Fig. 1: Flowchart of gene disruption by PCR-targeting



Protocol (see Flowchart Fig. 1)

1

Purification of the PCR template (resistance (-*oriT*) cassette)

Using whole plasmids as templates for the PCR can result in a high proportion of antibiotic-resistant transformants without gene disruption. This is caused by traces of CCC DNA that compete with the linear PCR fragment and result in the occurrence of false positive transformants. Using gel-purified disruption cassettes as templates prevents the occurrence of false positives.

1. Digest ~ 10 µg plasmid DNA (see Table 1) with 50 U *EcoRI* (Roche) and 50 U *HindIII* (Roche) in 1 X buffer B (Roche) in a 100 µl reaction.
 - A 2938 bp vector fragment and a fragment 14 bp larger than the size of the cassette given in Table 1 should be generated.
2. Run the digest on a 20 x 20 x 0.25 cm (100 ml) 1% TAE (1x) agarose gel at 5V/cm for 2 - 3 h in 1x TAE buffer.
 - Longer runs exhaust the buffer capacity and destroy the gel unless the buffer is recycled.
3. Cut out the cassette band from the gel and purify using the Qiagen gel extraction kit. The purified fragment is stored in 10 mM Tris.HCl (pH 8) at a concentration of 100 ng / µl at -20°C.
4. Absence of plasmid DNA is tested by using 1µl (100 ng) of purified cassette DNA to transform highly competent *E. coli* DH5α cells (10⁸/µg). Plate on LB agar containing 100 µg/ml carbenicillin. If any transformants appear, repeat steps 2-4.

2**Design of long PCR primers**


For each gene disruption, two long PCR primers (58 nt and 59 nt) are required. Each has at the 5' end 39 nt matching the *S. coelicolor* sequence adjacent to the gene to be inactivated, and a 3' sequence (19 nt or 20 nt) matching the right or left end of the disruption cassette (all cassettes have the same "right" and "left" ends). The precise positioning of the 39 nt sequence as indicated in Fig. 2 is important for creating in-frame deletions by FLP recombinase-induced excision of the resistance marker (see section 7).

- The 5'- 39 nt sequence of the forward primer (upstream primer; Fig. 2) must be from the coding strand of the gene of interest and its 3' end must be in the correct reading frame with respect to the replaced gene. The 5'- 39 nt sequence of the reverse primer (downstream primer; Fig. 2) must be from the *complementary* strand.
- To prevent unwanted recombination, a BlastN search is performed comparing each 39 nt sequence with the "real cosmid" (sequences at the Sanger Centre Homepage in the folder ftp.sanger.ac.uk/pub/S_coelicolor/cosmid_inserts and on the CD in the folder /S_coelicolor/cosmid inserts). The perfect match should be found but no other matches >30 bp. If necessary, the 39 nt sequence is shifted in 3 nt steps until the above criteria are met.


Fig.2: Designing PCR primers for making an *in-frame* deletion (the example illustrates a complete deletion)

39 nt from sense strand ending in ATG or GTG start codon
If sequence with >30 matches occurs in cosmid clone,
move → 3n nt to maintain frame

39 nt from anti-sense strand ending in Stop codon
← move -3n nt if necessary



FLP recombinase (BT340)



81 bp scar

(20bp + 19bp priming sequence + 42bp FLP core recombination site (see Fig.3); no in frame STOP)

3

PCR amplification of the extended resistance cassette

All PCR amplifications are performed using the Expand high fidelity PCR system according to the manufacturer’s instructions (Roche). Reaction conditions:

- Primers (100 pmoles/μl) 0.5 μl each 50 pmoles each
- Template DNA (100 ng/μl) 0.5 μl 50 ng ≈ 0.06 pmoles
- Buffer (10x) 5 μl 1 x
- dNTPs (10 mM) 1 μl each 50 μM each
- DMSO (100 %) 2.5 μl 5%
- DNA polymerase (2.5 U/μl) 1 μl 2.5 Units
- Water 36 μl
- Total volume 50 μl

Cycle conditions:

- 1. Denaturation: 94°C, 2 min
 - 2. Denaturation: 94°C, 45 sec
 - 3. Primer annealing: **50°C**, 45 sec
 - 4. Extension: 72°C, 90 sec
 - 5. Denaturation: 94°C, 45 sec
 - 6. Primer annealing: **55°C**, 45 sec
 - 7. Extension: 72°C, 90 sec
 - 8. Final extension: 72°C, 5 min
- } 10 cycles
- } 15 cycles

5 μl of the PCR product is used for analysis by gel electrophoresis. The expected sizes are 78 bp larger than the sizes of the disruption cassettes listed in Table 1 (because of the 2 x 39 bp 5’-primer extensions). The remaining 45 μl of the PCR product is purified using the Qiagen PCR purification kit according to the manufacturer’s instructions. The PCR product is finally eluted from the columns with 12 μl of water (~200 ng/μl).

4**Introduction of *S. coelicolor* cosmid clone into *E. coli* BW25113/pIJ790 (λ RED recombination plasmid) by electroporation**

pIJ790 contains the resistance marker *cat* (chloramphenicol resistance) and a temperature sensitive origin of replication (requires 30°C for replication).

1. Grow *E. coli* BW25113/pIJ790 overnight at 30°C in 10 ml LB (Luria-Bertani medium; Sambrook *et al.*, 1998) containing chloramphenicol (25 µg/ml).
2. Inoculate 100 µl *E. coli* BW25113/pIJ790 from overnight culture in 10 ml SOB (Hanahan, 1983) containing 20 mM MgSO₄ (add 200 µl of 1M stock to 10 ml SOB) and chloramphenicol (25 µg/ml).
3. Grow for 3-4 h at 30°C shaking at 200 rpm to an OD₆₀₀ of ~ 0.4.
4. Recover the cells by centrifugation at 4000 rpm for 5 min at 4°C in a Sorvall GS3 rotor (or equivalent).
5. Decant medium and resuspend the pellet by gentle mixing in 10 ml ice-cold 10 % glycerol.
6. Centrifuge as above and resuspend pellet in 5 ml ice-cold 10 % glycerol, centrifuge and decant. Resuspend the cell pellet in the remaining ~ 100 µl 10 % glycerol.
7. Mix 50 µl cell suspension with ~ 100 ng (1-2 µl) of cosmid DNA. Carry out electroporation in a 0.2 cm ice-cold electroporation cuvette using a BioRad GenePulser II set to: 200 Ω, 25 µF and 2,5 kV. The expected time constant is 4.5 – 4.9 ms.
8. Immediately add 1 ml ice cold LB to shocked cells and incubate shaking for 1h at 30°C.
9. Spread onto LB agar containing carbenicillin (100 µg/ml), kanamycin (50 µg/ml) and chloramphenicol (25 µg/ml).
10. Incubate overnight at 30°C.
11. Transfer one isolated colony into 5 ml LB containing antibiotics as in (9) above.
12. Incubate overnight at 30°C. This culture will be used as a pre-culture for generating competent cells to be transformed with the extended resistance cassette.

5

PCR targeting of the *S. coelicolor* cosmid

E. coli BW25113/pIJ790 containing a *S. coelicolor* cosmid is electro-transformed with the extended resistance cassette. **The example described uses the apramycin – *oriT* disruption cassette from pIJ773.** Table 1 lists alternative cassettes and their resistance determinants.

1. Inoculate a 10 ml SOB (*without* MgSO₄) culture containing carbenicillin (100 µg/ml), kanamycin (50 µg/ml) and chloramphenicol (25 µg/ml) with 1% of the overnight culture of *E. coli* BW25113/pIJ790 and the *S. coelicolor* cosmid. Add 100 µl 1M L-arabinose stock solution (final concentration is 10 mM, induces *red* genes).
2. Grow for 3-4 h at 30°C shaking at 200 rpm to an OD₆₀₀ of ~ 0.4.
3. Recover the cells by centrifugation at 4000 rpm for 5 min at 4°C in a Sorvall GS3 rotor (or equivalent).
4. Decant medium and resuspend the pellet by gentle mixing in 10 ml ice-cold 10% glycerol.
5. Centrifuge as above and resuspend pellet in 5 ml ice-cold 10 % glycerol, centrifuge and decant. Resuspend the cell pellet in remaining ~ 100 µl 10 % glycerol.
6. Mix 50 µl cell suspension with ~ 100 ng (1-2 µl) of PCR product. Carry out electroporation in a 0.2 cm ice-cold electroporation cuvette using a BioRad GenePulser II set to: 200 Ω, 25 µF and 2,5 kV. The expected time constant is 4.5 – 4.9 ms.
7. Immediately add 1 ml ice cold LB to shocked cells and incubated shaking 1 h at 37°C (or 30°C if further gene disruptions will be made on the same cosmid; see below).
8. Spread onto LB agar containing carbenicillin (100 µg/ml), kanamycin (50 µg/ml) and apramycin (50 µg/ml). If no further gene disruptions will be made on this cosmid, incubate overnight at 37°C to promote the loss of pIJ790. (If further disruptions are planned propagate overnight at 30°C and include chloramphenicol (25 µg/ml) so that pIJ790 is retained).

- If no colonies are obtained after 16 h growth at 37°C, repeat the experiment starting with a 50 ml SOB culture instead of 10 ml culture for generating electrocompetent cells. Try to concentrate the cells as much as possible by removing all of the remaining 10% glycerol. Resuspend the cell pellet in 50 µl 10% glycerol and use for electroporation.
- After 12 – 16 h growth at 37°C different colony-sizes are observed. Cultivating for longer time results in an increased background of small colonies, which are false positives. It is important to note that at this stage wild-type and mutant cosmids exist within one cell. The transformation with a PCR product and its integration in the cosmid DNA by homologous recombination will not occur in all copies of the cosmid molecules in one cell. One copy of a cosmid containing the incoming resistance marker is sufficient for resistance to this antibiotic. Normally, the larger the size of a colony, the more copies of mutagenised cosmids are present. Inoculating a large colony in 5 ml LB liquid cultures containing carbenicillin (100 µg/ml), kanamycin (50 µg/ml) and apramycin (50 µg/ml) result in a growth at 37°C to a cell density (OD₆₀₀ ~ 0.1 – 0.3) within 3-4 h (*E. coli* BW25113 without pIJ790 grows very fast). After 6 h plasmid DNA can be isolated and tested by restriction analysis and/or PCR using the primers described below.
- PCR analysis with a primer pair (test primers) priming just ~ 100 bp outside the region affected by homologous recombination will generate the expected fragment after gene disruption, but will usually also generate the wild-type fragment, caused by remaining wild-type copies within the same transformant. These will be lost during the subsequent transformation step into the methylation-deficient *E. coli* host ET12567 containing the non-transmissible plasmid pUZ8002 (this is not a problem anyway because wild-type copies lack the *oriT*).
- Notes on viomycin selection: selecting for viomycin^R depends critically on the amount of salt in the medium; more viomycin is required at higher salt concentrations. For a clean selection of *E. coli* clones, use DNA agar or 2xYT agar containing 30 µg/ml viomycin (see Kieser *et al.*, 2000).

For multiple gene replacements, choose an *oriT*-containing disruption cassette for the first knock-out, and a cassette without *oriT* and different resistance markers for further gene disruptions.

The gene disruption is confirmed by restriction analysis and/or PCR. Cosmid DNA of transformants is isolated from a 6 h, 37°C, 5 ml LB culture containing carbenicillin (100 µg/ml), kanamycin (50 µg/ml) and apramycin (50 µg/ml). Alkaline lysis followed by phenol/chloroform extraction produces cosmid DNA suitable for restriction analysis.

Cosmid CCC DNA isolation

1. Resuspend the cell pellet from 1 ml culture by vortexing in 100 µl solution I (50 mM Tris/HCl, pH 8; 10 mM EDTA).
2. Immediately add 200 µl solution II (200 mM NaOH; 1% SDS) and mix by inverting the tubes 10x.
3. Immediately add 150 µl solution III (3 M potassium acetate, pH 5.5) and mix by inverting the tubes 5x.
4. Spin at full speed in a microcentrifuge for 5 min at room temperature.
5. Immediately extract supernatant with 400 µl phenol/chloroform, vortex 2 min and spin at full speed in a micro centrifuge for 5 min.
6. Transfer the upper phase and add 600 µl 2-propanol. Leave the tubes on ice for 10 min.
7. Spin as above and wash the pellet with 200 µl 70% ethanol.
8. Spin as above and leave the tube open for 5 min at room temperature to dry the pellet. Resuspend the pellet in 50 µl 10mM Tris/HCl (pH 8) and use 10 µl for restriction digest.
 - **Omitting the phenol/chloroform extraction step results in degradation of the cosmid DNA. Use of miniprep-columns without including a phenol/chloroform extraction is not recommended.**

Verification of positive transformants by PCR requires an additional pair of 18 – 20 nt test primers which anneal 100 – 200 bp upstream and downstream of the 39 bp recombination region. (These primers can also be used later to verify the FLP-mediated excision of the resistance cassette.)

• Primers (100 pmoles/ μ l)	0.2 μ l each	20 pmoles each
• Template DNA (~50 ng/ μ l)	1 μ l	50 ng
• Buffer (10x)	5 μ l	1 x
• dNTPs (10 mM)	1 μ l each	50 μ M each
• DMSO (100 %)	2.5 μ l	5%
• DNA polymerase (2.5 U/ μ l)	1 μ l	2.5 Units
• Water	36.1 μ l	
• Total volume	50 μ l	

Cycle conditions:

1. Denaturation:	94°C, 2 min	} 30 cycles
2. Denaturation:	94°C, 45 sec	
3. Primer annealing:	55°C, 45 sec	
4. Extension:	72°C, 90 sec	
5. Final extension:	72°C, 5 min	

5 μ l of the PCR product is used for gel electrophoresis.

6

Transfer of the mutant cosmids into *Streptomyces*

If the target *Streptomyces* for mutagenesis carries a methyl-sensing restriction system (as is the case for *S. coelicolor* and *S. avermitilis*), it is necessary to passage the cosmid containing an apramycin resistance-*oriT* cassette through a non-methylating *E. coli* host. To achieve this, it is introduced by transformation into the non-methylating *E. coli* ET12567 containing the RP4 derivative pUZ8002. The cosmid is then transferred to *Streptomyces* by intergeneric conjugation (see Table 2 for resistance markers). If the target *Streptomyces* for mutagenesis does **not** carry a methyl-sensing restriction system (as is the case for *S. lividans*), common *E. coli* strains such as DH5 α containing pUZ8002 can be used instead.

Description	Name	Replication	Carb ^R	Cml ^R	Kan ^R	Tet ^R
<i>S. coelicolor</i> cosmid clones	Supercos 1		Carb ^R		Kan ^R	
λ Red plasmid	pIJ790	t ^s		Cml ^R		
FLP recombinase plasmid	BT340	t ^s	Carb ^R	Cml ^R		
OriT ⁻ RP4 derivative	pUZ8002				Kan ^R	
OriT ⁺ RP4 derivative	pUB307				Kan ^R	
Non-methylating <i>E. coli</i>	ET12567			Cml ^R		Tet ^R

Table 2. Resistance markers of vectors, helper plasmids and strains (carbenicillin resistance (Carb^R), chloramphenicol resistance (Cml^R), kanamycin resistance (Kan^R), tetracycline resistance (Tet^R), temperature sensitive replicon (t^s)). See Table 1 for replacement cassettes.

1. Prepare competent cells of *E. coli* ET12567/pUZ8002 grown at 37°C in LB containing kanamycin (25 μ g/ml) and chloramphenicol (25 μ g/ml) to maintain selection for pUZ8002 and the *dam* mutation, respectively. (ET12567 has a doubling time > 30 min.)
 - **High competence is required when Dam-methylated plasmids are introduced into a *dam*⁻ strain.**
2. Transform competent cells with the *oriT*-containing cosmid clone, and select for the incoming plasmid only using apramycin (50 μ g/ml) and carbenicillin (100 μ g/ml).

3. Inoculate a colony into 10 ml LB containing apramycin (50 µg/ml), chloramphenicol (25 µg/ml) and kanamycin (50 µg/ml). Grow overnight at 37°C.
 - **Chloramphenicol^S or Kanamycin^S segregants arise frequently among transformants, so set up more than one culture. The kanamycin selection is probably ineffective because both the cosmid and pUZ8002 confer resistance (Table 2).**
4. Inoculate 100 µl overnight culture into 10 ml fresh LB plus antibiotics as above and grow for ~ 4 h at 37°C to an OD₆₀₀ of 0.4.
5. Wash the cells twice with 10 ml of LB to remove antibiotics that might inhibit *Streptomyces*, and resuspend in 1 ml of LB.
6. While washing the *E. coli* cells, for each conjugation add 10 µl (10⁸) *Streptomyces* spores to 500 µl 2 × YT broth. Heat shock at 50°C for 10 min, then allow to cool.
7. Mix 0.5 ml *E. coli* cell suspension and 0.5 ml heat-shocked spores and spin briefly. Pour off most of the supernatant, then resuspend the pellet in the c. 50 µl residual liquid.
8. Make a dilution series from 10⁻¹ to 10⁻⁴ each step in a total of 100 µl of water.
9. Plate out 100 µl of each dilution on MS agar + 10mM MgCl₂ (without antibiotics) and incubate at 30°C for 16-20 h.
10. Overlay the plate with 1 ml water containing 0.5 mg nalidixic acid (20 µl of 25 mg/ml stock; selectively kills *E. coli*) and 1.25 mg apramycin (25 µl of 50 mg/ml stock). Use a spreader to lightly distribute the antibiotic solution evenly. Continue incubation at 30°C.
11. Replica-plate each MS agar plate with single colonies onto DNA plates containing nalidixic acid (25 µg/ml) and apramycin (50 µg/ml) with and without kanamycin (50 µg/ml). Double cross-over exconjugants are kanamycin^S and apramycin^R. (DNA gives fast, non-sporulating growth.)
12. Kanamycin^S clones are picked from the DNA plates and streaked for single colonies on MS agar (promotes sporulation) containing nalidixic acid (25 µg/ml) and apramycin (50 µg/ml).
13. Confirm kanamycin sensitivity by replica-plating onto DNA plates containing nalidixic acid (25 µg/ml) with and without kanamycin (50 µg/ml).
14. Purified kanamycin sensitive strains are then verified by PCR and Southern blot analysis.

- Typically, ~ 10 % of the exconjugants are double cross-over recombinants. The frequency of double cross-overs depends on the length of the flanking regions of homologous DNA on the cosmid. If < 1 kb is left on one side of the disrupted gene, obtaining kanamycin^S double cross-over types directly on the conjugation plates may be difficult. It may be necessary to streak out several exconjugants for single colonies on MS agar without antibiotics. After 3-5 days growth replica-plate onto DNA with and without kanamycin.

Antibiotic	Stock mg/ml	µl for 1 ml overlay	Final conc. after flooding µg/ml	Concentration in	
				MS, DNA µg/ml	R2YE µg/ml
Apramycin	50	25	50	50	50
Kanamycin	50	100	200	50	200
Spectinomycin	200	25	200	400	400
Streptomycin	10	25	10	10	10
Viomycin	30	25	30	30	NA
Nalidixic acid	25 in 0.3 M NaOH	20	20	25	25

Table 3: Antibiotic concentrations for selection on *S. coelicolor* MS conjugation plates, DNA replica plates or R2YE protoplast regeneration plates (Note some small differences from Kieser *et al.*, 2000).

7

FLP-mediated excision of the disruption cassette

The disruption cassettes are flanked by FRT sites (FLP recognition targets). Expression of the FLP-recombinase in *E. coli* removes the central part of the disruption cassette, leaving behind a 81 bp “scar” sequence which, in the preferred reading frame (bold in Fig. 3), lacks stop codons.

I	P	G	I	R	R	P	A	V	R	S	S	Y	S	L	E	S	I	G	T	S	K	Q	L	Q	P	T
F	R	G	S	V	D	L	Q	F	E	V	P	I	L	*	K	V	*	E	L	R	S	S	S	S	L	
S	G	D	P	S	T	C	S	S	K	F	L	F	S	R	K	Y	R	N	F	E	A	A	P	A	Y	
<u>ATTCGGGGATCCGTCGACCT</u>																										
	10		20			30			40			50			60			70			80					
TAAGCCCCCTAGGCAGCTGGACGTCAAGCTTCAAGGATAAGAGATCTTTCATATCCTTGAAG																										
N	R	P	D	T	S	R	C	N	S	T	G	I	R	*	F	T	Y	S	S	R	L	L	E	L	R	C
G	P	I	R	R	G	A	T	R	L	E	*	E	R	S	L	I	P	V	E	F	C	S	W	G	V	
E	P	S	G	D	V	Q	L	E	F	N	R	N	E	L	F	Y	L	F	K	S	A	A	G	A	*	

Fig. 3: Sequence of the 81 bp “scar” sequence remaining after FLP-mediated excision of the disruption cassette. The translation of the preferred reading frame is printed bold. The 20 and 19 nt priming sites are underlined and printed in colour. (Fig. 2 explains the determination of the reading frame.)

★ indicate stop codons,

priming site (20 nt)

priming site (19 nt)

This allows the generation of (hopefully) non-polar, unmarked in-frame deletions and repeated use of the same resistance marker for making multiple knock-outs in the same cosmid or in the same strain. *E. coli* DH5 α cells containing the temperature sensitive FLP recombination plasmid BT340 (Datsenko and Wanner, 2000; can be obtained from the *E. coli* Genetic Stock Center: CGSC Strain# 7629) are transformed with the mutagenised cosmid DNA (obtained in section 5). BT340 contains ampicillin and chloramphenicol resistance determinants and is temperature sensitive for replication (replicates at 30°C). FLP synthesis and loss of the plasmid are induced at 42°C (Cherepanov and Wackernagel, 1995).

1. Grow *E. coli* DH5 α /BT340 overnight at 30°C in 10 ml LB containing chloramphenicol (25 μ g/ml).
 - **Transforming *E. coli* BW25113/cosmid::apramycin (mutagenised cosmid) with the plasmid BT340 is not recommended because the isolates after PCR targeting may still contain copies of undisrupted cosmid DNA (see page 10, second paragraph).**
2. Inoculate 100 μ l *E. coli* DH5 α /BT340 from overnight culture into 10 ml LB containing chloramphenicol (25 μ g/ml).
3. Grow for 3-4 h at 30°C shaking at 200 rpm to an OD₆₀₀ of \sim 0.4.
4. Recover the cells by centrifugation at 4000 rpm for 5 min at 4°C in a Sorvall GS3 rotor (or equivalent).
5. Decant medium and resuspend the pellet by gentle mixing in 10 ml ice-cold 10 % glycerol.
6. Centrifuge as above and resuspend pellet in 5 ml ice-cold 10 % glycerol, centrifuge and decant. Resuspend the cell pellet in remaining \sim 100 μ l 10% glycerol.
7. Mix 50 μ l cell suspension with \sim 100 ng (1-2 μ l) of mutagenised cosmid DNA. Carry out electroporation in a 0.2 cm ice-cold electroporation cuvette using a BioRad GenePulser II set to: 200 Ω , 25 μ F and 2,5 kV. The expected time constant is 4.5 – 4.9 ms.
8. Immediately add 1 ml ice cold LB to shocked cells and incubate shaking for 1 h at 30°C.
9. Spread onto LB agar containing apramycin (50 μ g/ml) and chloramphenicol (25 μ g/ml).
10. Incubate for 2 d at 30°C (*E. coli* DH5 α /BT340 grows slowly at 30°C).
11. A single colony is streaked on an LB agar plate without antibiotics for single colonies and grown overnight at 42°C to induce expression of the FLP recombinase followed by the loss of plasmid BT340.
12. Make two masterplates by streaking 20 – 30 single colonies with a toothpick first on a LB agar plate containing apramycin (50 μ g/ml) and then on a LB agar plate containing kanamycin (50 μ g/ml).

13. Grow the masterplates overnight at 37°C. Apramycin^S kanamycin^R clones indicate the successful loss of the resistance cassette and are further verified by restriction and PCR analysis.

- Typically, ~ 10 % of the single colonies after non-selective growth lose the incoming resistance marker and the plasmid BT340 simultaneously.
- Using the same test primers as in section 5 (annealing ~ 100 bp upstream and downstream of the 39 nt primer sequence) should produce a PCR product of ~ 300 bp (200 bp + 81 bp “scar”). PCR fragments can be sequenced using the amplification primers for verification.

Replacing resistance cassette inserts in *S. coelicolor* with the unmarked “scar” sequence

The chromosomal apramycin resistance cassette insert in *S. coelicolor* is replaced by the “scar” sequence. This is achieved by homologous recombination between the chromosome and the corresponding “scar cosmid” prepared in 7. The procedure differs from section 6 because the cosmid lacks *oriT*, and the desired product is antibiotic sensitive. Therefore, it is necessary to introduce the scar cosmid into *Streptomyces* by protoplast transformation, and then select for kanamycin resistant *Streptomyces* containing the entire scar cosmid integrated by a single crossover. Restreaking to kanamycin-free medium, followed by screening for concomitant loss of kanamycin resistance and apramycin resistance, then identifies the desired *Streptomyces* clones.

Preparation of *Streptomyces coelicolor* protoplasts

1. Add 25 ml YEME medium to a baffled flask. Add ~ 0.1 ml spore suspension and required growth factors. Incubate 36-40 h at 30°C in an orbital incubator shaker.
 - Cultures of *S. lividans* and *S. coelicolor* are ready for harvesting when they start to produce red pigment
2. Pour culture broth into a 20 ml screw cap bottle and spin in the bench centrifuge (~ 1000 x g, 10 min).
 - Before centrifugation, examine the culture for contamination by unicellular bacteria, usually indicated by turbidity: the *Streptomyces* mycelium sediments

quickly while unicellular contaminants remain suspended. In case of doubt, use the microscope.

3. Discard the supernatant carefully; the pellet is easily disturbed.
 - **If the mycelium does not pellet add 5 ml sterile water to reduce the density of the medium and centrifuge again.**
4. Resuspend pellet in 15 ml 10.3% sucrose and spin in bench centrifuge as above. Discard supernatant.
5. Repeat step 4.
 - **The mycelial pellet, without added liquid, can be stored frozen at -20°C**
6. Resuspend mycelium in 4 ml lysozyme solution (1 mg/ml P buffer, filter sterilised); incubate at 30°C , 15-60 min.
7. Draw in and out of a 5 ml pipette three times and incubate for a further 15 min.
 - **This helps to free protoplasts from the mycelium so that they will pass through the cotton wool filter used in step 9. At least with *S. lividans*, it is possible to obtain transformants with unfiltered material, but the washing (steps 9-10) is still needed to remove lysozyme.**
8. Add 5 ml P buffer. Repeat step 7.
9. Filter protoplasts through cotton wool (using a filter tube) and transfer to a plastic tube.
10. Sediment protoplasts gently by spinning in a bench centrifuge ($\sim 1000 \times g$, 7 min).
11. Discard supernatant and suspend protoplasts in 1 ml P buffer.
 - **At this and any other steps when pelleted protoplasts are to be resuspended, resuspend in the remaining drop of liquid by tapping the side of the tube repeatedly with a finger until the protoplasts are dispersed to form a creamy suspension, then add the suspending P buffer (otherwise the protoplast pellet is difficult to disperse). Avoid vortexing, which induces foaming and consequent lysis. To freeze the protoplasts for storage, place samples of the protoplast suspension in small plastic tubes, close them and place them in ice in a plastic beaker. Place the beaker at -70°C overnight. Free the frozen protoplasts in their tubes from the ice and store at -70°C . To thaw, shake the frozen tube under running warm water (i.e. freeze slowly, thaw quickly). To assess the proportion of non-protoplasted units in the suspension, samples can be diluted in parallel in P buffer and in dilute detergent ($\sim 0.01\%$ SDS) and plated on regeneration plates. Any colonies arising after dilution in detergent are likely to have arisen from non-protoplasted units.**

Rapid small-scale transformation of *Streptomyces coelicolor*

1. Dispense 50 μl samples of protoplasts ($\sim 10^{10}/\text{ml}$) into as many tubes as there are transformations.
 - We usually spin the protoplasts down immediately before the transformation experiment. This eliminates substances that may have leaked out of the protoplasts during storage and the contents of protoplasts which have lysed spontaneously (which may include nucleases).
2. Complete steps 2a-c for each transformation individually.
 - a. Add up to 5 μl DNA solution to protoplasts and mix immediately by tapping tube.
 - b. Add 200 μl 25% PEG 1000 in P buffer and mix by pipetting up and down four times (be careful not to contaminate the barrel of the pipette).
 - c. Spread protoplast suspension (100-200 μl) on two dried R2YE plates. Use P buffer to make dilutions if required.
 - 1 ml glass pipettes can be used instead of spreaders. The solution will spread to some extent by itself if the plates are left on a horizontal surface.
3. Incubate plates at 30°C. After 14-20 h, flood for kanamycin selection. Score for resistant colonies after 3 d.
4. Select single colony and streak non-selectively for single colonies on MS agar plates and grow 3-4 d at 30°C.
5. Replica-plate to DNA agar plates with apramycin (50 $\mu\text{g}/\text{ml}$) or kanamycin (50 $\mu\text{g}/\text{ml}$) to screen for apramycin^S and kanamycin^S transformants.

The λ RED recombination plasmid pIJ790




Fig. 4: The λ RED recombination plasmid pKD20 (*E. coli* Genetic Sock Center CGSC Strain# 7637; Datsenko and Wanner, 2000) was modified by replacing the ampicillin resistance gene *bla* by the chloramphenicol resistance gene *cat*, generating pIJ790. Unlike pKD20, pIJ790 can be selected in presence of Supercos1.

Legend:

araC: encodes arabinose activator

bet: encodes single strand DNA binding protein

cat: chloramphenicol resistance gene

exo: exonuclease gene, promotes recombination together with *bet*

gam: gene product inhibits the host exonuclease V

oriR101: origin of replication

P araBAD: L-arabinose inducible promoter

repA101ts: temperature-sensitive replication

The following digestions can be used to verify this plasmid.

*Bam*HI: 6084 bp; *Eco*RI: 2872 bp, 1703bp, 1509 bp;

*Nco*I: 3927 bp, 2157 bp; *Pst*I: 5873 bp, 247 bp

Template plasmids




Fig. 5: Template plasmid pIJ773 containing the apramycin resistance gene *aac(3)IV* (AC=X99313) and the *oriT* of plasmid RP4 (=RK2) (AC=L27758), flanked by FRT sites (FLP recognition targets, see Datsenko and Wanner, 2000). The disruption cassette was cloned into the *EcoRV* site of pBluescript KS (+) allowing its isolation as a 1382 bp *EcoRI/HindIII* fragment.

* indicates suitable restriction sites for verification of mutagenised cosmid DNA by restriction analysis (for example: *SstI* generates a 751 bp internal fragment within the disruption cassette).

LOCUS pIJ773 4334 bp DNA CIRCULAR SYN 02-AUG-2002
 DEFINITION Ligation of Apra-oriT disruption cassette into the EcoRV site of
 pBluescript SK(+)
 ACCESSION pIJ773
 KEYWORDS .
 SOURCE Unknown.
 ORGANISM Unknown
 Unclassified.
 REFERENCE 1 (bases 1 to 4334)
 AUTHORS Gust et al., 2003,
 JOURNAL *Proc. Natl. Acad. Sci. USA* 100(4), 1541-1546
 FEATURES Location/Qualifiers
 CDS 2..20
 /region="priming site 19 bp"
 /product="TGT AGG CTG GAG CTG CTT C"
 CDS 21..54
 /region="FRT"
 /product="natural FRT site"
 CDS complement (66..851)
 /gene="aac(3)IV"
 /product="apramycin resistance gene"
 CDS 1192..1301
 /region="oriT"
 /product="origin of transfer (RK2)"
 CDS 1310..1343
 /region="FRT"
 /product="natural FRT site"
 CDS complement (1352..1371)
 /region="priming site 20 bp"
 /product="ATT CCG GGG ATC CGT CGA CC"
 CDS 2200..3060
 /gene="amp"
 /product="b-lactamase"
 BASE COUNT 1023 a 1137 c 1126 g 1049 t
 ORIGIN
 1 gtgtaggctg gagctgcttc gaagttccta tactttctag agaataggaa cttcggaata
 61 ggaacttatg agctcagcca atcgactggc gagcggcctc gcattcttcg catcccgcct
 121 ctggcggatg caggaagatc aacggatctc ggcccagttg acccagggct gtcgcccaca
 181 tgtcgcggga gcggatcaac cgagcaaagg catgaccgac tggaccttcc ttctgaaggg
 241 tcttctcctt gagccacctg tccggccaagg caaagcgctc acagcagtgg tcattctcga
 301 gataatcgac gcgtaccaac ttgccatcct gaagaatggt gcagtgtctc ggcaccccat
 361 agggaacctt tgccatcaac tcggcaagat gcagcgtcgt gttggcatcg tgtcccacgc
 421 cgaggagaag tacctgcccc tcgagttcat ggacacgggc gaccgggctt gcaggcgagt
 481 gaggtgcaag gggcaatgga tcagagatga tctgctctgc ctgtggcccc gctgcccga
 541 aggcfaatgg atgggctctg cgctttacat ttggcaggcg ccagaatgtg tcagagacaa
 601 ctccaaggtc cgggtgtaacg ggcgacgtgg caggatcgaa cggctcgtcg tccagacctg
 661 accacgaggg catgacgagc gtccctcccc gaccacggcg agcaccgagg gcctcgtatca
 721 ctccaagtgg cccatctctg agggcccgga cgctacggaa ggagctgtgg accacgaca
 781 caccgcccgg ggtaaccccc aggttgagaa gctgaccgat gagctcggct tttcgccatt
 841 cgtattgcac gacattgcac tccaccgctg atgacatcag tcgatcatag cacgatcaac
 901 ggcactgttg caaatagtcg gtgggtgataa acttatcatc cccttttgcg gatggagctg
 961 ccaatgaacc cattcaagg cggcattttt cagcgtgaca tcattctgtg gggcgtacgc
 1021 tggtaactgca aatacggcat cagttaccgt gagctgcatt ttccgctgca taaccctgct
 1081 tcggggctcat tatagcgatt ttttcggtat atccatcctt tttcgcacga tatacaggat
 1141 tttgcccagg ggttcgtgta gactttcctt ggtgtatcca acggcgtcag cggggcagga
 1201 taggtgaagt agggcccacc gcgacggggt gttccttctt cactgtccct tattcgcaac
 1261 tggcgggtgct caacgggaat cctgctctgc gaggtggcg ggaacttcga agttcctata
 1321 ctttctagag aataggaact tcgaactgca ggtcgacgga tcccgggaat atcaagctta
 1381 tcgataccgt cgacctcgag ggggggcccc gtaccaaat gcgccatag tgagtcgat
 1441 tacaattcac tggccgtcgt tttacaacgt cgtgactggg aaaaccctgg cgttacccaa
 1501 cttaatcgcc ttgcagcaca tccccctttc gccagctggc gtaatagcga agaggcccgc
 1561 accgatcgcc cttcccaca gttgcgcagc ctgaatggcg aatggcgcga cgcgccctgt
 1621 agcggcgcag taagcgcggc ggggtgtggg gttacgcgca gcgtgaccgc tacacttgcc
 1681 agcgcctag cgcccgtcc tttcgctttc ttcccttctt ttctcgccac gttcgccgc
 1741 tttcccgtc aagctctaaa tcgggggctc cctttagggt tccgatttag tgctttacgg
 1801 cacctcgacc ccaaaaaact tgattagggg gatggttcac gtagtggggc atcgccctga
 1861 tagacggttt ttgcgccctt gacggttgag tccacgttct ttaatagtgg actcttgctc
 1921 caaactggaa caacactcaa cctatctcg gtctattctt ttgattata agggattttg
 1981 ccgatttcgg cctattgggt aaaaaatgag ctgatttaac aaaaatttaa cgcgaatatt
 2041 aacaaaatat taacgtttac aatttcccag gtggcacttt tcggggaaat gtgcgaggaa
 2101 ccctattttg tttatttttc taaatacatt caaatatgta tccgctcatg agacaataac
 2161 cctgataaat gcttcaataa tattgaaaaa ggaagagtat gagtattcaa catttcctg
 2221 tcgcccctat tccccttttt cgggcatttt gccttctctg ttttctctac ccagaaacgc
 2281 tggtgaaagt aaaagatgct gaagatcagt tgggtgcac agtgggttac atcgactgg
 2341 atctcaacag cggtaagatc cttgagagtt ttcgccccga agaactgttt ccaatgatga
 2401 gcacttttaa agttctgcta tgtggcggcg tattatcccg tattgacgcc gggcaagagc
 2461 aactcggctg ccgcatacac tattctcaga atgacttggt tgagtactca ccagtcacag

2521 aaaagcatct tacggatggc atgacagtaa gagaattatg cagtgtgccc ataaccatga
2581 gtgataaacac tgcggccaac ttacttctga caacgatcgg aggaccgaag gagctaaccg
2641 cttttttgca caacatgggg gatcatgtaa ctgccttga tcggtgggaa cgggagctga
2701 atgaagccat accaaaacgac gagcgtgaca ccacgatgcc tgtagcaatg gcaacaacgt
2761 tgcgcaaaact attaactggc gaactactta ctctagcttc ccggcaacaa ttaatagact
2821 ggatggaggc ggataaagtt gcaggaccac ttctgctgc gccccctccg gctggctggt
2881 ttattgtctga taaatctgga gccggtgagc gtgggtctcg cggatcatt gcagcactgg
2941 ggccagatgg taagccctcc cgtatcgtag ttatctacac gacggggagt caggcaacta
3001 tggatgaacg aatagacag atcgcgtgaga taggtgcctc actgattaag cattggtaac
3061 tgtcagacca agtttactca tatatacttt agattgattt aaaacttcat ttttaattta
3121 aaaggatcta ggtgaagatc ctttttgata atctcatgac caaaatccct taacgtgagt
3181 tttcgttcca ctgagcgtca gaccccgtag aaaagatcaa aggatcttct tgagatcctt
3241 tttttctgcg cgtaatctgc tgcttgcaaa caaaaaaac accgctacca gcggtggttt
3301 gtttgccgga tcaagagcta ccaactcttt ttccgaaggt aactggcttc agcagagcgc
3361 agataccaaa tactgtcctt ctagtgtagc cgtagttagg ccaccacttc aagaactctg
3421 tagcaccgcc tacatacctc gctctgctaa tctgttacc agtggctgct gccagtggcg
3481 ataagtcgtg tcttaccggg ttggactcaa gacgatagtt accggataag gcgcagcggg
3541 cgggctgaac ggggggttcg tgcacacagc ccagcttggg gcgaacgacc tacaccgaac
3601 tgagatacct acagcgtgag ctatgagaaa gcgccacgct tcccgaaggg agaaaaggcg
3661 acaggtatcc ggtaagcggc agggctcggaa caggagagcg cacgagggag cttccagggg
3721 gaaacgcctg gtatctttat agtcctgtcg ggtttcgcca cctctgactt gagcgtcgat
3781 ttttgtgatg ctcgtcaggg gggcgagacc tatggaaaaa cgccagcaac gcggcctttt
3841 tacggttcct ggcttttgc tggccttttg ctcacatggt ctttctgcg ttatccctg
3901 attctgtgga taaccgtatt accgcctttg agtgagctga taccgctcgc cgcagccgaa
3961 cgaccgagcg cagcagatca gtgagcaggg aagcgggaaga gcgcccaata cgcaaaccgc
4021 ctctccccgc gcgttggccg attcattaat gcagctggca cgacagggtt cccgactgga
4081 aagcgggcag tgagcgaac gcaattaatg tgagttagct cactcattag gcaccccagg
4141 ctttacactt tatgcttccg gctcgtatgt tgtgtggaat tgtgagcggg taacaatttc
4201 acacagggaa cagctatgac catgattacg ccaagctcgg aattaaccct cactaaaggg
4261 aacaaaagct ggagctccac cgcggtggcg gccgctctag aactagtgga tccccgggg
4321 tgcaggaatt cgat




Fig. 6: Template plasmid pIJ778 containing the streptomycin and spectinomycin resistance gene *aadA* (AC=M60473) and the *oriT* of plasmid RP4 (=RK2) (AC=L27758), flanked by FRT sites (FLP recognition targets, see Datsenko and Wanner, 2000). The disruption cassette was cloned into the *EcoRV* site of pBluescript KS (+) allowing its isolation as a 1425 bp *EcoRI/HindIII* fragment.

* indicates suitable restriction sites for verification of mutagenised cosmid DNA by restriction analysis (for example: *NaeI* generates a 782 bp internal fragment within the disruption cassette).

LOCUS pIJ778 4377 bp DNA CIRCULAR SYN 02-AUG-2002

DEFINITION Ligation of Spec,Strep-oriT disruption cassette into the EcoRV site of pBluescript SK(+)

ACCESSION pIJ778

KEYWORDS .

SOURCE Unknown.

ORGANISM Unknown
Unclassified.

REFERENCE 1 (bases 1 to 4377)

AUTHORS Gust et al., 2003

JOURNAL Proc. Natl. Acad. Sci. USA 100(4), 1541-1546

FEATURES Location/Qualifiers

CDS 699..717
/region="priming site 19 bp"
/product="TGT AG complement (718.
/region="FRT"
/product="natural FRT site"

CDS complement (798..1590)
/gene="aadA"
/product="spectinomycin,streptomycin adenylyltransferase gene"

CDS 1931..2040
/region="oriT"
/product="origin of transfer (RK2)"

CDS 2049..2082
/region="FRT"
/product="natural FRT site"

CDS complement (2091..2110)
/region="priming site 20 bp"
/product="A . CC"

CDS complement (2111..2120)
/gene="amp"
/product="b-lactamase"

BASE COUNT 1065 a 1120 c 1108 g 1084 t

ORIGIN

1 GGAAATTGTA AACGTTAATA TTTTGTAAA ATTCGCGTTA AATTTTTGTT AAATCAGCTC
61 ATTTTTTAAC CAATAGGCCG AAATCGGCAA AATCCCTTAT AAATCAAAG AATAGACCGA
121 GATAGGGTTG AGTGTGTTC CAGTTTGGAA CAAGAGTCCA CTATTAAGA ACGTGGACTC
181 CAACGTCAA GGCAGAAA CCGTCTATCA GGCAGATGGC CCACTACGTG AACCATCACC
241 CTAATCAAGT TTTTGGGGT CGAGGTGCCG TAAAGCACTA AATCGGAACC CTAAAGGGAG
301 CCCCATTGAC AGAGCTTGAC GGGGAAAGCC GGCGAACGTG GCGAGAAAGG AAGGGAAGAA
361 AGCGAAAGGA GCGGGCGCTA GGGCGTGGC AAGTGTAGCG GTCACGCTGC GCGTAACCAC
421 CACACCGCC GCGCTTAATG CGCCGCTACA GGGCGGTCG CGCCATTCGC CATTAGGCT
481 GCGCAACTGT TGGGAAGGGC GATCGGTGCG GGCCCTCTCG CTATTACGCC AGCTGGCGAA
541 AGGGGATGT GCTGCAAGC GATTAAGTTG GGTAACGCCA GGGTTTCCC AGTCACGACG
601 TTGTAACACG ACGGCCAGTG AATTGTAATA CGACTCACTA TAGGGCGAAT TGGGTACCGG
661 GCCCCCCTC GAGGTCGACG GTATCGATAA GCTTGATGTG TAGGCTGGAG CTGCTCGAA
721 GTTCTATAC TTTCTAGAGA ATAGGAACCT CGGAATAGGA ACTTCATGAG CTCAGCCAAT
781 CGACTGGCGA GCGCATCTT ATTTGCCGAC TACCTTGGTG ATCTCGCTT TCACGTAGTG
841 GACAAATCT TCCAACGAT CTGCGCGCGA GGCCAAGCGA TCTTCTTCTT GTCCAAGATA
901 AGCCTGTCTA GCTCAAGTA TGACGGGCTG ATACTGGGCC GGCAGGCGCT CCATTGCCCA
961 GTCGGCAGC ACATCCTTCG GCGGATTTT GCGGTTACT GCGCTGTACC AAATGCGGGA
1021 CAACGTAAG ACTACATTC GCTCATGCC AGCCAGTCG GCGGGCGAGT TCCATAGCGT
1081 TAAGTTTTCA TTTAGCGCCT CAAATAGATC CTGTTCAGGA ACCGGATCAA AGAGTTCTCT
1141 CGCCGCTGGA CCTACCAAG CAACGCTATG TTCTCTTCTT TTTGTACGCA AGATAGCCAG
1201 ATCAATGTCG ATCGTGGCTG GCTCGAAGAT ACCTGCAAGA ATGTCAATTC GCTGCCATTC
1261 TCCAAATTCG AGTTCGCGCT TAGCTGGATA ACGCCACGGA ATGATGTCGT CGTGACAAAC
1321 AATGGTGACT TCTACAGCGC GGAGAATCTC GCTCTCTCCA GGGGAAAGCC AAGTTTCCAA
1381 AAGTTCGTTG ATCAAAGCTC GCCGCTTGT TTCATCAAGC CTTACGGTCA CCGTAACCAG
1441 CAAATCAATA TCACTGTGTG GCTTCAGGCC GCCATCCACT GCGGAGCCGT ACAAATGTAC
1501 GGCCAGCAAC GTCGGTTCGA GATGGCGCTC GATGACGCCA ACTACCTCTG ATAGTTGAGT
1561 CGATACTTCG GCGATCACCG CTTCCCTCAT GACATTGCAC TCCACCGCTG ATGACATCAG
1621 TCGATCATAG CACGATCAAC GGCAGTGTG CAAATAGTCG GTGGTGATAA ACTTATCATC
1681 CCCTTTTGCT GATGGAGCTG CACATGAACC CATTCAAAGG CCGGCATTTT CAGCGTGACA
1741 TCATTCTGTG GCGGTACGC TGGTACTGCA AATACGGCAT CAGTTACCGT GAGCTGCATT
1801 TTCCGTGCA TAACCTGCT TCGGGTTCAT TATAGCGATT TTTTCGGTAT ATCCATCCTT
1861 TTTCGCACGA TATACAGGAT TTTGCCAAAG GGTTCGTGTA GACTTTCCTT GGTGTATCCA
1921 ACGGCGTCA CCGGGCAGGA TAGGTGAAGT AGGCCACCC CCGAGCGGGT GTTCTTCTT
1981 CACTGTCCCT TATTGCGACC TGGCGGTGCT CAACGGGAAT CCTGCTCTGC GAGGCTGGCG
2041 GAGCCCTCGA AGTTCCTATA CTTTCTAGAG AATAGGAAT TCGAACTGCA GGTTCGACGGA
2101 TCCCCGGAAT ATCGAATTC TGCAGCCCGG GGGATCCACT AGTTCCTAGAG CCGCCGCCAC
2161 CGCGGTGGAG CTCAGCTTT TGTCCCTTT AGTGAGGGTT AATTCAGGAG TTGGCGTAAT
2221 CATGGTCATA GCTGTTTCCT GTGTGAAATT GTTATCCGCT CACAATFCCA CACAACATAC
2281 GAGCCGGAAG CATAAAGTGT AAAGCCTGGG GTGCCTAATG AGTGAGCTAA CTCACATFCA
2341 TTGCGTTGCG CTCAGTCCC GCTTTCAGT CCGGAAACCT GTCGTGCCAG CTGCATTAAT
2401 GAATCGGCCA ACGCGCGGG AGAGCGGTT TCGTATTGG GCGCTCTCC GCTTCTCCG

2461 TCACTGACTC GCTGCGCTCG GTCGTTCCGGC TCGGGCGAGC GGTATCAGCT CACTCAAAGG
2521 CGGTAATACG GTTATCCACA GAATCAGGGG ATAACGCAGG AAAGAACATG TGAGCAAAAG
2581 GCCAGCAAAA GGCAGGAAC CGTAAAAAGG CCGCGTTGCT GGCGTTTTTC CATAGGCTCC
2641 GCCCCCTGA CGAGCATCAC AAAAATCGAC GCTCAAGTCA GAGGTGGCGA AACCCGACAG
2701 GACTATAAAG ATACCAGGCG TTCCCCCTG GAAGCTCCCT CGTGCGCTCT CCTGTTCCGA
2761 CCCTGCCGCT TACCGGATAC CTGTCCGCCT TTCTCCCTTC GGAAGCGTG GCGCTTTCTC
2821 ATAGCTCACG CTGTAGGTAT CTCAGTTCGG GTAGGTCGT TCGCTCCAAG CTGGGCTGTG
2881 TGCACGAACC CCCGTTTCCAG CCCGACCGCT GCGCCTTATC CGGTAACTAT CGTCTTGAGT
2941 CCAACCCGCT AAGACACGAC TTATCGCCAC TGGCAGCAGC CACTGGTAAC AGGATTAGCA
3001 GAGCGAGGTA TGTAGGCGGT GCTACAGAGT TCTTGAAGTG GTGGCCTAAC TACGGCTACA
3061 CTAGAAGGAC AGTATTTGGT ATCTGCGCTC TGCTGAAGCC AGTTACCTTC GGAAAAAGAG
3121 TTGGTAGCTC TTGATCCGGC AAACAAACCA CCGCTGGTAG CCGTGGTTTT TTTGTTTGCA
3181 AGCAGCAGAT TACGCGCAGA AAAAAGGAT CTCAAGAAGA TCCTTTGATC TTTTCTACGG
3241 GGTCTGACGC TCAGTGGAAC GAAAACTCAC GTTAAGGGAT TTTGGTCATG AGATTATCAA
3301 AAAGGATCTT CACCTAGATC CTTTTAAATT AAAAATGAAG TTTTAAATCA ATCTAAAGTA
3361 TATATGAGTA AACTTGGTCT GACAGTTACC AATGCTTAAT CAGTGAGGCA CCTATCTCAG
3421 CGATCTGTCT ATTTCGTTCA TCCATAGTTG CCTGACTCCC CGTCGTGTAG ATAACTACGA
3481 TACGGGAGGG CTTACCATCT GGCCCCAGTG CTGCAATGAT ACCGCGAGAC CCACGCTCAC
3541 CGGCTCCAGA TTTATCAGCA ATAAACCAGC CAGCCGGAAG GGCCGAGCGC AGAAGTGGTC
3601 CTGCAACTTT ATCCGCCTCC ATCCAGTCTA TTAATTGTTG CCGGGAAGCT AGAGTAAGTA
3661 GTTCGCCAGT TAATAGTTTG CGCAACGTTG TTGCCATTGC TACAGGCATC GTGGTGTCCAC
3721 GCTCGTCGTT TGGTATGGCT TCATTACGCT CCGGTTCCCA ACGATCAAGG CGAGTTACAT
3781 GATCCCCCAT GTTGTGCAAA AAAGCGGTTA GTCCTTCGG TCCTCCGATC GTTGTACAGAA
3841 GTAAGTTGGC CGCAGTGTTA TCACTCATGG TTATGGCAGC ACTGCATAAT TCTCTTACTG
3901 TCATGCCATC CGTAAGATGC TTTTCTGTGA CTGGTGAGTA CTCAACCAAG TCATTCTGAG
3961 AATAGTGTAT GCGGCGACCG AGTTGCTCTT GCCCGGCGTC AATACGGGAT AATACCGCGC
4021 CACATAGCAG AACTTTAAAA GTGCTCATCA TTGGAAAACG TTCTTCGGGG CGAAAACTCT
4081 CAAGGATCTT ACCGCTGTTG AGATCCAGTT CGATGTAACC CACTCGTGCA CCCAACTGAT
4141 CTTACGCATC TTTTACTTTC ACCAGCGTTT CTGGGTGAGC AAAAACAGGA AGGCAAAATG
4201 CCGCAAAAAA GGAATAAAGG GCGACACGGA AATGTTGAAT ACTCATACTC TTCTTTTTTC
4261 AATATTATTG AAGCATTAT CAGGGTTATT GTCTCATGAG CGGATACATA TTTGAATGTA
4321 TTTAGAAAAA TAAACAAATA GGGGTTCCGC GCACATTTC CCGAAAAGTG CCACCTG

//




Fig. 7: Template plasmid pIJ780 containing the viomycin resistance gene *vph* (AC=X99314) and the *oriT* of plasmid RP4 (=RK2) (AC=L27758), flanked by FRT sites (FLP recognition targets, see Datsenko and Wanner, 2000). The disruption cassette was cloned into the *EcoRV* site of pBluescript KS (+) allowing its isolation as a 1497 bp *EcoRI/HindIII* fragment.

* indicates suitable restriction sites for verification of mutagenised cosmid DNA by restriction analysis (for example: *HincII* generates a 1137 bp internal fragment within the disruption cassette).

LOCUS pIJ780 4449 bp DNA CIRCULAR SYN 02-AUG-2002
 DEFINITION Ligation of Vio-oriT disruption cassette into the EcoRV site of
 pBluescript SK(+)
 KEYWORDS .
 SOURCE Unknown.
 ORGANISM Unknown
 Unclassified.
 REFERENCE 1 (bases 1 to 4449)
 AUTHORS Gust et al., 2003
 JOURNAL Proc. Natl. Acad. Sci. USA 100(4), 1541-1546
 FEATURES Location/Qualifiers
 CDS 2..20
 /region="priming site 19 bp"
 /product="TGT AGG CTG GAG CTG CTT C"
 CDS 21..54
 /region="FRT"
 /product="natural FRT site"
 CDS complement (102..925)
 /gene="vph"
 /product="viomycin resistance gene"
 CDS 1306..1415
 /region="oriT"
 /product="origin of transfer (RK2)"
 CDS 1424..1457
 /region="FRT"
 /product="natural FRT site"
 CDS complement (1466..1485)
 /region="priming site 20 bp"
 /product="ATT CCG GGG ATC CGT CGA CC"
 CDS 2314..3174
 /gene="amp"
 /product="b-lactamase"
 BASE COUNT 977 a 1241 c 1194 g 1037 t
 ORIGIN
 1 GTGTAGGCTG GAGCTGCTTC GAAGTTCCTA TACTTTCTAG AGAATAGGAA CTTCGGAATA
 61 GGAACTTCAT GAGCTCAGCC AATCGACTGG CGAGCGGCAT CCTACCGGTA GCCGCTGAGG
 121 CCGTCGGCGA GTTCTCCTC GTCGCCGTCG CGCTGCGCGT AGAGGGCCTG CTGGAGTGGC
 181 AAGGTGCCTC GGATCGCCGA GATCCGCTCG GCCGTTCCGT TGTCGGCCCA GCCGCCGAGC
 241 GCGAGCACTC GGCCAGCAG TTCCTCGCCG TAGCTCGCCC CGATGGCGGC CAGGTCCTCA
 301 GCCGGTGC C GATGCCGAC CTCGTCCCAG TCGACGACGC CGCTCATGCG CGGCACCTCCG
 361 TCCACCGTCT CCCACAGGAC GTTCTCGCCG CCGAGGTCAC CGTGGACCAC CGCGGAGGTG
 421 AGATGGGGCA GGGCGTCGAG CGCGCGGAGC TCGCGCTCGG CACGCTCCCG GCCCGCTCG
 481 GACATCAGCG GGAACAGTTC GGTACGCACC CCCGTGGCGA ACTCCTGCCA CTCGTTGCGG
 541 GGAGCCTCCG GCAGCGCGGC GCGCACCTTC TCCTCGTCGC CCGCCCGCCG GAGCCCGGAC
 601 AGCAGGGTCG CGTACTGTCG GCGCAGCGCC TCCGCCACCT CCGGGCTGGT GAGCACATCG
 661 TCCTCCAACG GTGCTCCGGG AATGCGGCTC AGCACCAAGT ACGGCGGCTC GTCCCTGCC
 721 TGGGCGCCGC CCTCGGACAG CGGCTGCGGC GTGCGAAACC CGAGGTCGAT CCCGGCAAGA
 781 GCGCGCAGGA CGTCCGCCCT GCCGGGCGA CGGTGCGCGG CCGCCCGGGT GCGGGCGAAG
 841 CAGACCACCC GGTGCGATCC GATCACCACA TGGTGGAACT GCCCTCGTG GACGGCGAGT
 901 CCGCCACCGG TGTCGCCCGG CAGGAGCCGG CTCAGCAGAT CGCGGTGCGT CTCAAATGATT
 961 CTCATGACAT TGCACTCCAC CGCTGATGAC ATCAGTCGAT CATAGCACGA TCAACGGCAC
 1021 TGTTGCAAAT AGTCGGTGGT GATAAACTTA TCATCCCTTT TTGCTGATGG AGCTGCACAT
 1081 GAACCCATTC AAAGGCCGGC ATTTTCAGCG TGACATCATT CTGTGGGCCG TACGCTGGTA
 1141 CTCGCAAATC GGCATCAGTT ACCGTGAGCT GCATTTTCCG CTGCATAACC CTGCTCGGG
 1201 GTCATTATAG CGATTTTTTC GGTATATCCA TCCTTTTTTC CACGATATAC AGGATTTTTG
 1261 CAAAGGGTTC GTGTAGACTT TCCTTGGTGT ATCCAACGGC GTCAGCCGGG CAGGATAGGT
 1321 GAAGTAGGCC CACCCGCGAG CGGGTGTTC TTCTTCACTG TCCCTTATTC GCACCTGGCG
 1381 GTGCTCAACG GGAATCCTGC TCTGCGAGGC TGCGCGGAAC TTCGAAGTTC CTATACTTTC
 1441 TAGAGAATAG GAACTTCGAA CTGCAGGTCG ACGGATCCCC GGAATATCAA GCTTATCGAT
 1501 ACCGTCGACC TCGAGGGGGG GCCCGGTACC CAATTGCCCC TATAGTGAGT CGTATACAAA
 1561 TTCACTGGCC GTCGTTTTAC AACGTCGTGA CTGGGAAAAC CCTGGCGTTA CCCAACTTAA
 1621 TCGCCTTGCA GCACATCCCC CTTTCGCCAG CTGGCGTAAT AGCGAAGAGG CCCGCACCGA
 1681 TCGCCCTTCC CAACAGTTGC GCAGCCTGAA TGGCGAATGG CGCGACGCGC CCTGTAGCGG
 1741 CGCATAAAGC GCGCGGGGTG TGGTGGTTAC GCGCAGCGTG ACCGCTACAC TTGCCAGCGC
 1801 CCTAGCGCCC GCTCCTTTCG CTTTCTTCCC TTCTTTTCTC GCCACGTTTC CCGGCTTTC
 1861 CCGTCAAAGCT CTAATCCGGG GGCTCCCTTT AGGGTTCGGA TTTAGTGCTT TACGCAACCT
 1921 CGACCCAAA AAACCTGATT AGGGTGATGG TTCACGTAGT GGGCCATCGC CCTGATAGAC
 1981 GGTTTTTCCG CCTTTGACGT TGGAGTCCAC GTTCTTTAAT AGTGGACTCT TGTTCCAAAC
 2041 TGGAACAACA CTCACCCCTA TCTCGGTCTA TTCTTTTGAT TTATAAGGGA TTTTGGCCGT
 2101 TTCGGCCTAT TGGTTAAAAA ATGAGCTGAT TTAACAAAAA TTTAACGCGA ATTTTAAACA
 2161 AATATTAACG TTTACAATTT CCCAGGTGGC ACTTTTCGGG GAAATGTGCG CGGAACCCCT
 2221 ATTTGTTTTT TTTTCTAAAT ACATTCAAAT ATGTATCCGC TCATGAGACA ATAACCTGTA
 2281 TAAATGCTTC AATAATATTG AAAAAGGAAG AGTATGAGTA TTCAACATTT CCGTGTGCGG
 2341 CTTATTCCCT TTTTTCGCGC ATTTTGCCTT CCGTGTTTTG CTCACCCAGA AACGCTGGTC
 2401 AAAGTAAAAG ATGCTGAAGA TCAGTTGGGT GCACGAGTGG GTTACATCGA ACTGGATCTC

2461 AACAGCGGTA AGATCCTTGA GAGTTTTTCGC CCCGAAGAAC GTTTTCCAAT GATGAGCACT
2521 TTTAAAGTTC TGCTATGTGG CGCGGTATTA TCCCGTATTG ACGCCGGGCA AGAGCAAATC
2581 GGTCCGCCGA TACACTATTC TCAGAATGAC TTGGTTGAGT ACTCACCAGT CACAGAAAAG
2641 CATCTTACGG ATGGCATGAC AGTAAGAGAA TTATGCAGTG CTGCCATAAC CATGAGTGAT
2701 AACACTGCGG CCAACTTACT TCTGACAACG ATCGGAGGAC CGAAGGAGCT AACCGCTTTT
2761 TTGCACAACA TGGGGGATCA TGTAACTCGC CTTGATCGTT GGGAAACCGGA GCTGAATGAA
2821 GCCATACCAA ACGACGAGCG TGACACCACG ATGCCTGTAG CAATGGCAAC AACGTTGCGC
2881 AAACATTTAA CTGGCGAACT ACTTACTCTA GCTTCCCAGC AACAAATTAAT AGACTGGATG
2941 GAGGCGGATA AAGTTGCAGG ACCACTTCTG CGTCGGGCC TTCCGGCTGG CTGGTTTATT
3001 GCTGATAAAT CTGGAGCCGG TGAGCGTGGG TCTCGCGGTA TCATTGCAGC ACTGGGGCCA
3061 GATGGTAAGC CCTCCCCTAT CGTAGTTATC TACACGACGG GGAGTCAGGC AACTATGGAT
3121 GAACGAAATA GACAGATCGC TGAGATAGGT GCCTCACTGA TTAAGCATTG GTAAGTGTCA
3181 GACCAAGTTT ACTCATATAT ACTTTAGATT GATTTAAAAC TTCATTTTTA ATTTAAAAGG
3241 ATCTAGGTGA AGATCCTTTT TGATAATCTC ATGACCAAAA TCCCTTAAAC TGAGTTTTTCG
3301 TTCACCTGAG CGTCAGACCC GTAGAAAAAG ATCAAAGGAT CTTCTTGAGA TCCTTTTTTT
3361 TTGCGCGTAA TCTGCTGCTT GCAAAACAAA AAACCACCGC TACCAGCGGT GGTGTTTTTG
3421 CCGATCAAG AGTACCAAC TCTTTTTCCG AAGGTAACG GCTTCAGCAG AGCGCAGATA
3481 CCAAATACTG TCCTTCTAGT GTAGCCGTAG TTAGGCCACC ACTTCAAGAA CTCTGTAGCA
3541 CCGCCTACAT ACCTCGCTCT GCTAATCCTG TTACCAGTGG CTGCTGCCAG TGGCGATAAG
3601 TCGTGCTTTA CCGGGTTGGA CTCAAGACGA TAGTTACCGG ATAAGGCGCA GCGGTCGGGC
3661 TGAACGGGGG GTTCGTGCAC ACAGCCCAGC TTGGAGCGAA CGACCTACAC CGAACTGAGA
3721 TACCTACAGC GTGAGCTATG AGAAAGCGCC ACGCTTCCCG AAGGGAGAAA GGCGGACAGG
3781 TATCCGGTAA GCGGCAGGGT CGGAACAGGA GAGCGCACGA GGGAGCTTCC AGGGGAAAC
3841 GCCTGGTATC TTTATAGTCC TGTCGGGTTT CGCCACCTCT GACTTGAGCG TCGATTTTTG
3901 TGATGCTCGT CAGGGGGCGG GAGCCTATGG AAAAACGCCA GCAACGCGGC CTTTTTACGG
3961 TTCTGGCCT TTTGCTGGCC TTTTGCTCAC ATGTTCTTTC CTGCGTTATC CCCTGATTCT
4021 GTGATAACC GTATTACCGC CTTTGAGTGA GCTGATAACG CTCGCCGAG CCGAACGACC
4081 GAGCGCAGCG AGTCAGTGAG CGAGGAAGCG GAAGAGCGCC CAATACGCAA ACCGCCTCTC
4141 CCCGCGGTT GCGCGATTCA TTAATGCAGC TGGCACGACA GGTGTTCCGA CTGGAAGCG
4201 GGCAGTGAGC GCAACGCAAT TAATGTGAGT TAGCTCACTC ATTAGGCACC CCAGGCTTTA
4261 CACTTTATGC TTCCGGCTCG TATGTTGTGT GGAATTGTGA GCGGATAACA ATTTACACA
4321 GGAAACAGCT ATGACCATGA TTACGCCAAG CTCGGAATTA ACCCTCACTA AAGGGAACAA
4381 AAGCTGGAGC TCCACCGCGG TGGCGGCCG TCTAGAATA GTGGATCCCC CGGGCTGCAG
4441 GAATTCGAT

//

Strains:

Strain designation: BW25113/pIJ790

Plasmid: pIJ790 [*oriR101*], [*repA101(ts)*], *araBp-gam-be-exo*

Chromosome: (Δ (*araD-araB*)567, Δ *lacZ4787(::rrnB-4)*, *lacIp-4000(lacI^Q)*, λ , *rpoS369(Am)*, *rph-1*, Δ (*rhaD-rhaB*)568, *hsdR514*; Datsenko and Wanner, 2000);

this strain should be grown on rich medium containing chloramphenicol (25 μ g/ml) at 30°C.

Strain designation: DH5 α /pIJ773

Plasmid: pBluescript KS (+), *aac(3)IV*, *oriT* (RK2), FRT sites

chromosomal marker : see Stratagene,

this strain should be grown on rich medium containing carbenicillin (100 μ g/ml) and apramycin (50 μ g/ml) at 37°C.

Strain designation: DH5 α /pIJ778

Plasmid: pBluescript KS (+), *aadA*, *oriT* (RK2), FRT sites

chromosomal marker : see Stratagene,

this strain should be grown on rich medium containing carbenicillin (100 μ g/ml), streptomycin (50 μ g/ml) and spectinomycin (50 μ g/ml) at 37°C.

Strain designation: DH5 α /pIJ780

Plasmid: pBluescript KS (+), *vph*, *oriT* (RK2), FRT sites

chromosomal marker : see Stratagene,

this strain should be grown on DNA (Difco Nutrient Agar), 2XYT broth or DNB (Difco Nutrient broth) containing carbenicillin (100 μ g/ml) and viomycin (30 μ g/ml) at 37°C.

Strain: DH5 α /BT340 (Datsenko and Wanner, 2000) can be ordered from:

E. coli Genetic Stock Center

830 Kline Biology Tower

MCD Biology Department

266 Whitney Ave. Box 208103

Yale University

New Haven, CT 06520-81033

Tel. (0044) (203) 432-9997 Fax –6161

email: berlyn@cgsc.biology.yale.edu

webpage: <http://cgsc.biology.yale.edu>

Strains DH5 α /pIJ779 and DH5 α /pIJ781 will be send separately on request,

please contact:

Plant Bioscience Limited, Dr Karin Schofield

Technology Acquisition & Technology Manager Norwich Research Park, Colney,

Norwich, NR4 7UH, UK Email: karin@plantbioscience.com

Tel: +44 (0) 1603 456500; Fax: +44 (0) 1603 456552

www.plantbioscience.com

Media: for more detailed information see:

Kieser T, Bibb MJ, Buttner MJ, Chater KF and Hopwood DA (2000)

Practical Streptomyces Genetics, John Innes Foundation, Norwich Research Park,

Colney, Norwich NR4 7UH, England

Agar Media

Note: MS medium (Mannitol Soya flour Medium) is also known as SFM medium (Kieser *et al.*, 2000).

R2 Medium (Okanishi *et al.*, 1974; Hopwood and Wright, 1978)

Make up the following solution:

Sucrose	103 g
K ₂ SO ₄	0.25 g
MgCl ₂ .6H ₂ O	10.12 g
Glucose	10 g
Difco Casaminoacids	0.1 g
Distilled water	800 ml

Pour 80 ml of the solution into 250 ml Erlenmeyer flasks each containing 2.2 g Difco Bacto agar. Close the flasks and autoclave. At time of use, remelt the medium and add to each flask the following autoclaved solutions in the order listed:

KH ₂ PO ₄ (0.5%)	1 ml
CaCl ₂ .2H ₂ O (3.68%)	8 ml
L-proline (20%)	1.5 ml
TES buffer (5.73%, adjusted to pH7.2)	10 ml
¹ Trace element solution	0.2 ml
NaOH (1N) (unsterilised is OK)	0.5 ml
Required growth factors for auxotrophs	0.75 ml each

(For stock solutions, see Kieser *et al.*, 2000)

¹Trace element solution (l⁻¹):

ZnCl ₂	40 mg
FeCl ₃ .6H ₂ O	200 mg
CuCl ₂ .2H ₂ O	10 mg
MnCl ₂ .4H ₂ O	10 mg
Na ₂ B ₄ O ₇ .10H ₂ O	10 mg
(NH ₄) ₆ Mo ₇ O ₂₄ .4H ₂ O	10 mg

R2YE Medium

Make up in the same way as R2, but add 5 ml of Difco yeast extract (10%) to each 100 ml flask at time of use

R5 Medium

This is an alternative way of making R2YE. It obviously saves time at the bench. For most purposes it seems satisfactory, but there is some folklore suggesting that R2YE made by the traditional process gives more reproducible results.

Make up the following solution:

Sucrose	103 g
K ₂ SO ₄	0.25 g
MgCl ₂ .6H ₂ O	10.12 g
Glucose	10 g
Difco Casaminoacids	0.1 g
¹ Trace element solution	2 ml
Difco yeast extract	5 g
TES buffer	5.73 g
Distilled water to	1000 ml

¹ The same as in R2 and R2YE.

Pour 100 ml of the solution into 250 ml Erlenmeyer flasks each containing 2.2 g Difco Bacto agar. Close the flasks and autoclave. At time of use, re-melt the medium and add to each flask in the order listed:

KH ₂ PO ₄ (0.5%)	1 ml
CaCl ₂ .2H ₂ O (5M)	0.4 ml
L-proline (20%)	1.5 ml
NaOH (1N)	0.7 ml
Any required growth factors for auxotrophs (see Kieser <i>et al.</i> , 2000)	0.75 ml each

Mannitol Soya flour Medium (MS) Hobbs *et al.* (1989). Sometimes referred to as “SFM”.

Agar	20 g
Mannitol	20 g
¹ Soya flour	20 g
Tap water	1000 ml

¹Use soya flour from a health food shop or supermarket, not the expensive material from a laboratory supplier.

Dissolve the mannitol in the water and pour 200 ml into 250 ml Erlenmeyer flasks each containing 2 g agar and 2 g soya flour. Close the flasks and autoclave *twice* (115 °C , 15 min), with gentle shaking between the two runs.

Difco nutrient agar (DNA)

Place 4.6 g Difco Nutrient Agar in each 250 ml Erlenmeyer flask and add 200 ml distilled water. Close the flasks and autoclave.

L agar

Agar	10 g
Difco Bacto tryptone	10 g
NaCl	5 g
Glucose	1 g
Distilled water	up to 1000 ml

Dissolve the ingredients, except agar, in the distilled water and pour 200 ml into 250 ml Erlenmeyer flasks each containing 2 g agar. Close the flasks and autoclave.

Liquid media

Yeast extract-malt extract medium (YEME)

Difco yeast extract	3 g
Difco Bacto-peptone	5 g
Oxoid malt extract	3 g
Glucose	10 g
Sucrose	340 g (34% final)
Distilled water	up to 1000 ml

After autoclaving, add:

MgCl ₂ .6H ₂ O (2.5M)	2 ml/litre (5mM final)
---	------------------------

For preparing protoplasts, also add:

¹ Glycine (20%)	25 ml/litre (0.5% final)
----------------------------	--------------------------

¹Different *Streptomyces* strains may need different concentrations of glycine; 0.5% is best for *S. lividans* and *S. coelicolor*.

Tryptone soya broth (TSB)

Oxoid Tryptone Soya Broth powder (CM129)	30 g
Distilled water	1000 ml

Difco nutrient broth (DNB)

Difco Nutrient Broth powder	8 g
Distilled water	1000 ml

L (Lennox) broth (LB)

Difco Bacto tryptone	10 g
Difco yeast extract	5 g
NaCl	5 g
Glucose	1 g
Distilled water	1000 ml

2 X YT medium

Difco Bacto tryptone	16 g
Difco Bacto yeast extract	10 g
NaCl	5 g
Water	1000 ml

Buffers

P (protoplast) buffer

Make up the following basal solution:

Sucrose	103 g
K ₂ SO ₄	0.25 g
MgCl ₂ .6H ₂ O	2.02 g
¹ Trace element solution	2 ml
Distilled water to	800 ml

¹ The same as in R2 and R2YE.

Dispense in 80 ml aliquots and autoclave. Before use, add to each flask in order:

KH ₂ PO ₄ (0.5%)	1 ml
CaCl ₂ .2H ₂ O (3.68%)	10 ml
TES buffer (5.73%, adjusted to pH7.2)	10 ml

T (transformation) buffer

Mix the following sterile solutions:

Sucrose (10.3%)	25 ml
Distilled water	75 ml
¹ Trace element solution	0.2 ml
K ₂ SO ₄ (2.5%)	1 ml

¹ The same as in R2 and R2YE.

To 9.3 ml of the above solution add:

CaCl ₂ (5M)	0.2 ml
¹ Tris-maleic acid buffer	0.5 ml

¹Make up a 1M solution of Tris and adjust to pH8.0 by adding maleic acid. For use, add 3 parts by volume of the above solution to 1 part by weight of PEG 1000, previously sterilised by autoclaving.

L (lysis) buffer (Thompson *et al.*, 1982)

Mix the following sterile solutions:

Sucrose (10.3%)	100 ml
TES buffer (5.73%, adjusted to pH7.2)	10 ml
K ₂ SO ₄ (2.5%)	1 ml
¹ Trace element solution	0.2 ml
KH ₂ PO ₄ (0.5%)	1 ml
MgCl ₂ ·6H ₂ O (2.5M)	0.1 ml
CaCl ₂ (0.25M)	1 ml

¹ The same as in R2 and R2YE.

This stock solution keeps indefinitely. Just before use dissolve lysozyme in a sample of the solution at a concentration of 1 mg/ml and sterilise by filtration.

The CD-ROM includes:

1. Manual as Microsoft Word document and Adobe Acrobat document (X:/Protocol/)
2. Sequences of template plasmids as plain sequence files and as EMBL files (X:/Template plasmids/)
3. Sequences of the *S. coelicolor* real cosmids (X:/S_coelicolor/cosmid inserts)

Literature

- Bennett, P.M., Grinsted, J. and Richmond, M.H. (1977) *Mol.Gen.Genet.*, 154, 205-211**
- Cherepanov, P.P and Wackernagel, W. (1995) *Gene*, 158, 9-14**
- Datsenko, K.A. and Wanner, B.L. (2000) *Proc. Natl. Acad. Sci. USA* 97, 6640-6645**
- Flett, F., Mersinias, V. and Smith, C.P. (1997) *FEMS Microbiol.Lett.*, 155, 223-229**
- Hobbs G, Frazer CM, Gardner DCJ, Cullum JA, Oliver SG (1989) *Appl.Microbiol.Biotechnol.*, 31, 272-277.**
- Hodgson DA (1982) *J.Gen.Microbiol.*,128, 2417-2430.**
- Hopwood DA (1967) *Bacteriol.Rev.*, 31, 373-403.**
- Hopwood DA, Bibb MJ, Chater KF, Kieser T, Bruton CJ, Kieser HM, Lydiat DJ, Smith CJ, Ward JM, and Schrepf H. (1985) Norwich: The John Innes Foundation**
- Kieser, T., Bibb, M.J., Buttner, M.J., Chater, K.F. and Hopwood, D.A. (2000) *Practical Streptomyces Genetics*: John Innes Foundation, Norwich Research Park, Colney, Norwich NR4 7UH, UK**
- MacNeil, D.J., Gewain, K.M., Ruby, C.L., Dezeny, G., Gibbons, P.H. and MacNeil, T. (1992) *Gene*, 111, 61-68**
- Matsushima ,P., Broughton, M.C., Turner, J.R. and Baltz, R.H. (1994) *Gene*, 146, 39-45**
- Redenbach, M., Kieser, H.M., Denapaite, D., Eichner, A., Cullum, J., Kinashi, H. and Hopwood, D.A. (1996) *Mol. Microbiol.*, 21, 77-96**
- Sambrook, J., Fritsch, E.F. and Maniatis, T. (1989) Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.**
- Strauch E, Takano E, Baylis H, Bibb MJ (1991) *Mol.Microbiol.*, 5, 289-298.**